
## **COOPERATIVE LIGAND BINDING: OXYGEN BINDING TO HEMOGLOBIN**

## Lab Activity

- 1. Create 3 mL of a 1:200 dilution of bovine blood with 50 mM phosphate buffer (pH 7.4) in a small glass test tube.
- 2. Blank the spectrophotometer with 50 mM phosphate buffer (pH 7.4).
- 3. Discard the buffer in the cuvette, and transfer the diluted blood into the cuvette.
- 4. Add a small magnetic stir bar.
- 5. Connect the cuvette to the variable-pressure apparatus:



- 6. Record atmospheric gas pressure in mm Hg.
- 7. Record the current room temperature.
- 8. Through the two-way valve, pull the vacuum to 35 mm Hg, being careful not to let the contents of the cuvette boil. Close the two-way valve.
- 9. Position the cuvette on a magnetic stir plate and gently stir the contents of the cuvette for 30 minutes.
- 10. Position the cuvette in the spectrophotometer and read the absorbance between 520 and 590 nm.
- 11. Increase the pressure inside the variable-pressure cuvette apparatus. To add air in a controlled manner, fill a 20 mL syringe with a volume of 10 mL air and attach the syringe to the two-way value. Open the value, allowing the air to enter the system. Close the value, and remove the syringe. Record the new gas pressure.
- 12. Position the cuvette on a magnetic stir plate and gently stir the contents of the cuvette for 2 minutes.
- 13. Position the cuvette in the spectrophotometer and read the absorbance between 520 and 590 nm.

14. Repeat steps 11-13 until the gas pressure inside the system has returned to atmospheric pressure.

## Data Analysis

1. Determine the  $pO_2$  present in the system from the absolute pressure of air within the apparatus ( $P_{observed}$ ) for each observed pressure reading:

Temperature (°C)	pH <sub>2</sub> O (mm Hg)
18	15.5
19	16.5
20	17.5
21	18.7
22	19.8

$$pO_2 = 0.21 \left( P_{observed} - pH_2O \right)$$

- 2. Make a plot of  $A_{541 nm}$  versus pO<sub>2</sub>. Is hemoglobin behaving non-cooperatively or cooperatively?
- 3. Use non-linear regression with Microsoft Excel to solve for the baseline, amplitude, p50, and h parameters of the Hill model for your data at  $A_{541 nm}$ :

$$A_{\text{observed}} = \text{baseline} + \text{amplitude} \left( \frac{\left( pO_2 \right)^h}{\left( P_{50} \right)^h + \left( pO_2 \right)^h} \right)$$

- 4. What are the standard deviations associated with each parameter?
- 5. Add the calculated  $A_{541 nm}$  values to the plot created for Question #2 to compare the experimental data to the best-fit curve.