The standard Gibbs free energy change for converting 3-phosphoglycerate to 2-phosphoglycerate is 1.12 kcal mole⁻¹. Within a fed cardiac myocyte, the Gibbs free energy change for converting 3phosphoglycerate to 2-phosphoglycerate is -0.14 kcal mole⁻¹.

- 1. At equilibrium, is/3-phosphoglycerate or 2-phosphoglycerate molecule is more prevalent?
- Within a fed cardiac myocyte, is the conversion of 3-phosphoglycerate to 2-phosphoglycerate unfavorable or unfavorable?
- 3. Within a fed cardiac myocyte, is the conversion of 3-phosphoglycerate to 2-phosphoglycerate spontaneous or nonspontaneous?
- 4. For the condition within a fed cardiac myocyte, complete the following expression: